

ILLUSTRATED PARTS LIST

for the


Kodak DirectView CLASSIC/ELITE CR SYSTEMS

Service Code: 5258, 5259



Important

Qualified service personnel must repair this equipment.



H219_0001GC

PLEASE NOTE The information contained herein is based on the experience and knowledge relating to the subject matter gained by Carestream Health, Inc. prior to publication.

No patent license is granted by this information.

Carestream Health, Inc. reserves the right to change this information without notice, and makes no warranty, express or implied, with respect to this information. Carestream Health shall not be liable for any loss or damage, including consequential or special damages, resulting from any use of this information, even if loss or damage is caused by Carestream Health's negligence or other fault.



This equipment includes parts and assemblies sensitive to damage from electrostatic discharge. Use caution to prevent damage during all service procedures.

Table of Contents

Description	Page
Illustrated Parts List	3
SLOW SCAN AY	3
OPTICAL AY	5
Electrical Components	7
PMT COLLECTOR AY	9
PANELS and COVERS	11
Internal PANELS	13
SUPPORT PLATE AY - Version 1	15
ERASE AY	17
ANTI-REVERSE AY	19
CASSETTE HANDLING AY - Left	21
CASSETTE HANDLING AY - Right	23
BAR CODE LABELS	25
CABLES	27
Special TOOLS	29
FLOOR STAND	31
WALL STAND	33
DESKTOP ACCESSORIES	35
Numerical Index	37
Alphabetical Index	45
Publication History	53

Section 1: Illustrated Parts List

Figure 1 SLOW SCAN AY





Figure 1 SLOW SCAN AY

Item	Part No.	Description	Quantity	Notes
1	SP7H6068	SLIDE AY - slow scan, complete	1	
2	SP6H9899	SENSOR AY - plate present	1	
3	SP6H3359	EXTRACTION BAR AY	1	
4	SP943600	SCREW - <i>Torx</i> , m4 x 18 mm	2	Includes: • Y-HOOK MOTOR AY • INTERNAL BAR CODE READER
5	SP7H6238	MOTOR AY - extraction bar	1	
6	---	SLIDE AY - slow scan	1	
7	---	MOTOR GUARD - slow scan	1	
8	---	GROUND STRAP	1	
9	SP748360	SCREW - <i>Torx</i> , m4 x 12 mm	8	
10	---	BOLT - shoulder	1	
11	---	FLAT WASHER - 0.375 x 0.170 x 0.024 in.	2	
12	---	SLIDE - encoder AY	1	
13	5E1927	COUPLING - rotor	1	
14	3E1275	ROTOR AY	1	
15	9F4123	BOARD - slow scan, coil	1	
16	3E1274	GASKET - thermal	1	
17	---	BEARING - needle	1	
18	SP6H4562	PLATE - slider, extraction bar	1	
19	SP6H4656	SPRING - compression, extraction bar	1	
20	SP6H6768	CAP - extraction bar	1	
21	SP6H4565	PLATE - hook, extraction bar	1	
22	SP6H7815	SPRING - extension, extraction bar	1	
23	SP6H9388	STATIC BRUSH - extraction bar	1	
24	1F2412	BAR CODE READER - internal	1	
25	SP6H4575	MOTOR AY - y-hook	1	

03DEC07
IP5258-1
Page
5 of 53

Figure 2 OPTICAL AY



H219_1002DCA
H219_1002DC

Figure 2 OPTICAL AY

Item	Part No.	Description	Quantity	Notes
1	SP7H7435	MINI FAST SCAN - complete	1	
2	1F5679	FILTER - laser, ND laser safety	1	
3	---	FILTER BRACKET - laser	1	
4	---	THUMBSCREW - m3	1	
5	SP748361	SCREW - <i>Sems</i> , hex flange head	1	
6	1F8430	GALVO AY	1	
7	SP7H7434	MINI FAST SCAN	1	Does not include the GALVO AY.
8	---	COVER - optics ay	1	

03DEC07
IP5258-1
Page
7 of 53

Figure 3 Electrical Components





Figure 3 Electrical Components

Item	Part No.	Description	Quantity	Notes
1	SP6H3377	BOARD - mini mcb	1	
2	SP7H5533	BOARD - IEB	1	
3	---	RACK AY - power supply	1	
4	---	CONNECTOR	1	
5	SP7H3043	TRANSFORMER - isolation	1	
6	SP7H2777	HARNESS - power supply	1	
7	3E1359	RELAY	1	
8	---	STANDOFF	1	
9	---	COVER - relay	1	
10	---	PLUG	1	
11	SP7H3310	POWER SUPPLY	1	
12	---	BRACKET - power supply	1	
13	477506	CARRIER - fuse	1	
14	597341	FUSE - 8A, 250 V, slow blow	1	
15	943366	HOLDER - fuse	1	
16	SP7H6000	CIRCUIT BREAKER - power supply	1	
17	SP6H8080	SWITCH - main power	1	
18	SP7H2860	SUPPORT - transformer	1	

03DEC07
IP5258-1
Page
9 of 53

Figure 4 PMT COLLECTOR AY




H219_1003DCA
H219_1003DC

Figure 4 PMT COLLECTOR AY

Item	Part No.	Description	Quantity	Notes
1	---	COVER - pmt	1	
2	SP7H6142	COLLECTOR AY - complete	1	
3	---	COLLECTOR AY - 2 PMT	1	
4	5H7918	PMT AY	2	
5	1F4044	FOAM - pmt	2	
6	SP7E7744	SCREW - <i>Sems</i> , hex, m4 x 10 mm	2	
7	---	HOUSING - PMT	1	
8	SP7H3855	BOARD - PMT, 1.0 mz	1	
9	---	MOUNT - collector, rear	1	
10	SP943600	SCREW - <i>Torx</i> , m4 x 18 mm	2	
11	6H1776	FILTER - fan	1	
12	SP6H9897	FAN AY - collector	1	
13	---	MOUNT - flat cable	5	
14	SP7H5420	FOLD MIRROR AY	1	
15	---	MOUNT - fold mirror	1	
16	---	FOLD MIRROR	1	
17	---	MOUNT - collector, front	1	
18	SP6H8426	BAFFLE - afterglow		

03DEC07
IP5258-1
Page
11 of 53

Figure 5 PANELS and COVERS




H219_1000DCA
H219_1000DC

Figure 5 PANELS and COVERS

Item	Part No.	Description	Quantity	Notes
1	SP6H5471	COVER - top	1	
2	SP7H5993	LIGHT SEAL AY - left	1	
3	SP7H5992	LIGHT SEAL AY - right	1	
4	SP6H5474	NOSE PIECE - top cover	1	
5	SP6H4591	COVER - left	1	
6	SP6H4957	COVER - front	1	
7	SP6H4590	COVER - right	1	
8	SP7H4386	SCREW - captive	8	
9	SP6H6527	CASSETTE GUIDE - right	1	
10	SP6H6528	CASSETTE GUIDE - left	1	
11	SPSP2B4081	SCREW - <i>Torx</i> , m4 x 8 mm	4	
12	SP7H3856	CONTROL DISPLAY AY	1	

03DEC07
IP5258-1
Page
13 of 53

Figure 6 Internal PANELS


03DEC07
IP5258-1
Page
14 of 53

Figure 6 Internal PANELS


Item	Part No.	Description	Quantity	Notes
1	---	SEAL - light-power switch	2	
2	---	EMI PANEL - left	1	
3	---	EMI PANEL - rear	1	
4	---	SCREW - <i>Torx</i> , m4 x 18 mm	1	
5	---	EMI PANEL - right	1	
6	SP6H8005	GASKET - EMI, erase, long	2	
7	SP6H8007	GASKET - EMI, erase, short	2	
8	SP6H4975	EMI PANEL AY - front	1	
9	3H3775	LATCH AY - erase	1	

03DEC07
IP5258-1
Page
15 of 53

Figure 7 SUPPORT PLATE AY - Version 1



1




H219_1005DCA
H219_1005DC

03DEC07
IP5258-1
Page
16 of 53

Figure 7 PLATE SUPPORT AY - Version 1

Item	Part No.	Description	Quantity	Notes
1	SP7H6652	SUPPORT AY - plate	1	
2	---	BALL BEARING - flanged, 0.25 x 0.37 in.	2	
3	---	WASHER - curved, 0.49 x 0.25 x 0.011 in.	1	
4	---	STANDOFF - clinch, 10-32	1	
5	---	RIVET - 0.123 x 7/32 in.	1	
6	---	SHOULDER SCREW - 10-32 x 0.25 x 1.00 in.	1	
7	---	WASHER - 0.62 x 0.255 x 0.062 in.	2	
8	---	PIVOT BRACKET - plate support	1	
9	---	SKID - arm, plate support	1	
10	---	ARM - plate support	1	

03DEC07
IP5258-1
Page
17 of 53

Figure 8 ERASE AY

H219_1006DCA
H219_1006DC

Figure 8 ERASE AY

Item	Part No.	Description	Quantity	Notes
1	SP7H2685	ERASE AY - complete	1	
2	7E9471	ERASE LAMP	10	
3	SP5F7653	SUPPORT - lamp	1	
4	SP7H2968	BAFFLE AY - erase ay	1	
5	---	LABEL - erase ay	1	
6	SP182590	WIRE TIE	1	
7	365159	MOUNT - wire tie	1	
8	SP9B9445	SCREW - <i>Torx</i> m2.5 x 10 mm	1	
9	1C0554	TRANSFORMER - erase ay	1	
10	943134	SCREW - <i>Torx</i> m4 x 8 mm	1	
11	---	HARNESS - ballasts and sockets	1	
12	---	HARNESS - plug and sensor	1	
13	SP748360	SCREW - <i>Torx</i> , m4 x 12 mm	1	
14	---	FRAME - erase ay	1	
15	914047	LOCK NUT - hex	1	
16	---	WASHER - lock	1	
17	324259	SCREW - <i>Phillips</i> , M3 x 0.5 x 5 mm	1	
18	---	WASHER - lock, m3, internal tooth	1	
19	731068	WASHER - m3	1	

03DEC07
IP5258-1
Page
19 of 53

Figure 9 ANTI-REVERSE AY




Figure 9 ANTI-REVERSE AY

Item	Part No.	Description	Quantity	Notes
1	SP6H7946	ANTI-REVERSE AY - plate support	1	
2	---	POST PIVOT - eccentric	2	
3	---	ARM BOOSTER - plate support	1	
4	---	EXTENSION SPRING - arm booster	1	
5	---	BRACKET - disc, friction	1	
6	---	LEAF SPRING - anti-reverse	1	
7	---	SCREW - <i>Torx</i> , m4 x 6 mm	2	
8	---	SHOULDER SCREW - M4	1	
9	---	SHOULDER SCREW - M4 x 5.0 x 25 mm	2	
10	---	COMPRESSION DISC - outer	1	
11	---	WASHER - 0.62 x 0.255 x 0.062 in.	1	
12	---	CLUTCH - 22 x 8 x 9 mm	1	
13	---	PUSHER ARM - plate support	1	
14	SP6H9808	EXTENSION SPRING - clutch	1	
15	---	DISC AY - plate support	1	
16	---	COMPRESSION DISK - inner	1	
17	1F6393	SPRING - extension, anti-reverse ay	1	

03DEC07
IP5258-1
Page
21 of 53

Figure 10 CASSETTE HANDLING AY - Left





Figure 10 CASSETTE HANDLING AY - Left

Item	Part No.	Description	Quantity	Notes
1	SP7H3363	CASSETTE HANDLING AY - complete	1	
2	SP7H4840	HARNESS - mini-mcb to cassette handling ay	1	
3	SP7H3840	SPRING - light seal ay, left	1	
4	851123	SCREW - 10 - 32 x 0.375 in.	2	
5	852961	LOCK WASHER - light	2	
6	SP325608	DRIVE BELT - cassette handling ay	1	
7	SP7H6765	SPRING - extension, clamp	1	
8	---	SPRING PIN - 0.094 x 0.812 in., steel	1	
9	---	PULLEY - drive motor, 20T	1	
10	SP9F0341	MOTOR - cassette handling ay, DC	1	
11	SP7H2391	DRIVE MOTOR AY - cassette handling ay	1	
12	SP6H5260	PLUSH AY - left	1	
13	699579	SENSOR - emitter	1	
14	699577	SENSOR - detector	1	
15	TL3E1279	SWITCH - optical	1	
16	1F6158	PLUSH - light seal, front	1	
17	SP6H5836	SHUTTER AY - front	1	
18	SP6H6544	RAMP - entrance	1	
19	SP6H6063	DRIVE AY	1	

03DEC07
IP5258-1
Page
23 of 53

Figure 11 CASSETTE HANDLING AY - Right




H219_1009HCA
H219_1009HC

Figure 11 CASSETTE HANDLING AY - Right

Item	Part No.	Description	Quantity	Notes
1	SP7H3991	PLUSH - 750, right	1	
2	1F5817	PLUSH - 750, left	1	
3	6E5215	LIGHT SEAL - PLUSH, rear	5	
4	---	BRACKET - plush, rear	1	
5	SP6H6814	REAR CLIP AY	1	
6	K1J7774	SENSOR - GAP, wide	1	
7	---	SENSOR BRACKET - rear	1	
8	SP6H5190	LIGHT SHUTTER AY - right	1	
9	SP7H4100	PLUSH - clamp bar	2	
10	SP6H6423	SPRING - torsion, right plush	1	
11	SP6H6991	WEAR BAR - clamp, left	1	
12	1F7768	WEAR BAR - clamp, right	1	

03DEC07
IP5258-1
Page
25 of 53

Figure 12 BAR CODE LABELS



H187_0825HCA
H187_0825HC

Figure 12 BAR CODE LABELS

Item	Part No.	Description	Quantity	Notes
1	7C6678	LABEL - 35 x 43 cm, GP, 136 speed	1	
	1F2423	LABEL - 35 x 43 cm, LLI, 136 speed	2	
	5H0496	LABEL - 35 x 43 cm, GP+, 136 speed	2	
	7C6677	LABEL - 24 x 30 cm, GP, 136 speed	2	
	7C6680	LABEL - 24 x 30 cm, HR, 68 speed	2	
	9F3240	LABEL - 24 x 30 cm, EHR-M, 80 speed	2	
	1F2450	LABEL - 24 x 30 cm, HR-M, 68 speed	2	
	7C6676	LABEL - 18 x 24 cm, GP, 136 speed	2	
	7C6679	LABEL - 18 x 24 cm, HR, 68 speed	2	
	1F2451	LABEL - 18 x 24 cm, HR-M, 68 speed	2	
	9F3239	LABEL - 18 x 24 cm, EHR, 80 speed	2	
	9B8733	LABEL - 35 x 35 cm, GP, 136 speed	2	
	5H0495	LABEL - 35 x 35 cm, GP+, 136 speed	2	
	1F2424	LABEL - 15 x 30 cm, dental, 136 speed	1	
	6E5893	LABEL - tube side	2	
2 3	7E2130	LABEL - storage phosphor warning	2	
	1C7904	LATCH AY - cassette	1	For all CASSETTES.
	955887	SPRING - latch ay, cassette	1	For all CASSETTES.

03DEC07
IP5258-1
Page
27 of 53

Figure 13 CABLES





Figure 13 CABLES

Item	Part No.	Description	Quantity	Notes
1	SP7H4840	HARNESS - mini-mcb to cassette handling ay	1	
2	SP6H8030	CABLE - PMT	1	
3	SP6H8543	CABLE - laser control	1	
4	SP6H8029	CABLE - galvo motor	1	
5	---	HARNESS - erase ay	1	
	SP6H8594	CABLE - ieb board to mini-mcb board	1	Not visible in graphic.
	SP6H8593	CABLE - ethernet	1	Not visible in graphic.
	SP7H2640	CABLE - VGA	1	Not visible in graphic.

03DEC07
IP5258-1
Page
29 of 53

Figure 14 Special TOOLS




H219_1011HCA
H219_1011HC

03DEC07
IP5258-1
Page
30 of 53

Figure 14 Special TOOLS

Item	Part No.	Description	Quantity	Notes
	TL7H4608	ADJUSTMENT TOOL KIT		Includes all TOOLS.
1	---	ADJUSTMENT TOOL - shutter AY	1	
2	---	ADJUSTMENT TOOL - clamp alignment	1	
3	---	ADJUSTMENT TOOL - motor AY	1	
4	---	ADJUSTMENT TOOL - sensor S3	1	
5	---	ADJUSTMENT TOOL - bar code label	1	
6	---	ADJUSTMENT TOOL - fold mirror	1	
7	---	ADJUSTMENT TOOL - static brush	1	

03DEC07
IP5258-1
Page
31 of 53

Figure 15 FLOOR STAND


03DEC07
IP5258-1
Page
32 of 53

Figure 15 FLOOR STAND

Item	Part No.	Description	Quantity	Notes
1	---	STAND AY - monitor	1	
2	---	BIN - cassette storage	1	
3	914047	LOCK NUT - hex	7	
4	---	SHELF - monitor	1	
5	SP7H4879	TRAY - keyboard, floor stand	1	
6	---	LEVELING FOOT - 1.50 x 2.24 in.	3	
7	---	BASE - floor stand	1	
8	---	BRACKET - floor stand	1	
9	SP7H5101	LINER AY - bin, floor stand	1	

03DEC07
IP5258-1
Page
33 of 53

Figure 16 WALL STAND



H219_1014DCA
H219_1014DC

Figure 16 WALL STAND

Item	Part No.	Description	Quantity	Notes
1	SP7H5083	TRAY - keyboard, wall stand	1	
2	---	MONITOR	1	
3	SP7H5107	MONITOR ARM - wall stand	1	
4	---	COMPUTER BRACKET	1	
5	---	FRONT COMPUTER BRACKET	1	
6	---	THUMBSCREW - computer bracket	1	
7	---	MOUNTING BRACKET - wall stand	2	
8	---	WASHER	2	
9	---	SCREW	2	
10	---	WALL STAND BRACKET	1	
11	SP7H5089	ARM - bar code	1	
12	---	BRACKET - bar code reader	1	
13	SP7H2956	BOOT - bar code reader, wall stand	1	
14	SP7H5399	BAR CODE READER - handheld	1	
	SP7H5400	USB CABLE - bar code reader	1	

03DEC07
IP5258-1
Page
35 of 53

Figure 17 DESKTOP ACCESSORIES




Figure 17 DESKTOP ACCESSORIES

Item	Part No.	Description	Quantity	Notes
1	---	KEYBOARD - standard	1	Included with COMPUTER.
2	SP7H5087	KEYBOARD - mini	1	
3	SP6H7880	LCD MONITOR - 19 in., touchscreen	1	
	SP7H2639	POWER SUPPLY - 19 in. monitor	1	
	SP7H2641	CABLE - 19 in. monitor, USB	1	
	SP7H2642	CABLE - 19 in. monitor, audio	1	
4	SP7H3964	STAND - 19 in. monitor	1	Includes KEYBOARD and MOUSE, does not include software.
	SP6H7883	LCD MONITOR - 17 in.	1	
5	SP6H9581	COMPUTER - model 8810	1	
6	SP7H5401	STAND - handsfree bar code reader	1	
7	SP7H5399	BAR CODE READER - handheld	1	
	SP7H5400	USB CABLE - bar code reader	1	
	SP7H5403	CUP - bar code reader, desktop	1	
8	SP3B1966	POWER CORD - europe	1	
9	955815	POWER CORD - US&C hospital grade	1	
10	SP4B0826	POWER CORD - UK	1	
11	SP5E5632	POWER CORD - japan	1	
12	SP5E6714	POWER CORD - china	1	

Section 2: Numerical Index

Numerical Index

Part No.	Description	Figure No.
---	SLIDE AY - slow scan	1
---	MOTOR GUARD - slow scan	1
---	GROUND STRAP	1
---	BOLT - shoulder	1
---	FLAT WASHER - 0.375 x 0.170 x 0.024 in.	1
---	SLIDE - encoder AY	1
---	BEARING - needle	1
---	FILTER BRACKET - laser	2
---	THUMBSCREW - m3	2
---	COVER - optics ay	2
---	RACK AY - power supply	3
---	CONNECTOR	3
---	STANDOFF	3
---	COVER - relay	3
---	PLUG	3
---	BRACKET - power supply	3
---	COVER - pmt	4
---	COLLECTOR AY - 2 PMT	4
---	HOUSING - PMT	4
---	MOUNT - collector, rear	4
---	MOUNT - flat cable	4
---	MOUNT - fold mirror	4
---	FOLD MIRROR	4
---	MOUNT - collector, front	4
---	SEAL - light-power switch	6
---	EMI PANEL - left	6
---	EMI PANEL - rear	6
---	SCREW - Torx, m4 x 18 mm	6
---	EMI PANEL - right	6
---	BALL BEARING - flanged, 0.25 x 0.37 in.	7
---	WASHER - curved, 0.49 x 0.25 x 0.011 in.	7

Numerical Index

Part No.	Description	Figure No.
---	STANDOFF - clinch, 10-32	7
---	RIVET - 0.123 x 7/32 in.	7
---	SHOULDER SCREW - 10-32 x 0.25 x 1.00 in.	7
---	WASHER - 0.62 x 0.255 x 0.062 in.	7
---	PIVOT BRACKET - plate support	7
---	SKID - arm, plate support	7
---	ARM - plate support	7
---	LABEL - erase ay	8
---	HARNESS - ballasts and sockets	8
---	HARNESS - plug and sensor	8
---	FRAME - erase ay	8
---	WASHER - lock	8
---	WASHER - lock, m3, internal tooth	8
---	POST PIVOT - eccentric	9
---	ARM BOOSTER - plate support	9
---	EXTENSION SPRING - arm booster	9
---	BRACKET - disc, friction	9
---	LEAF SPRING - anti-reverse	9
---	SCREW - <i>Torx</i> , m4 x 6 mm	9
---	SHOULDER SCREW - M4	9
---	SHOULDER SCREW - M4 x 5.0 x 25 mm	9
---	COMPRESSION DISC - outer	9
---	WASHER - 0.62 x 0.255 x 0.062 in.	9
---	CLUTCH - 22 x 8 x 9 mm	9
---	PUSHER ARM - plate support	9
---	DISC AY - plate support	9
---	COMPRESSION DISK - inner	9
---	SPRING PIN - 0.094 x 0.812 in., steel	10
---	PULLEY - drive motor, 20T	10
---	BRACKET - plush, rear	11
---	SENSOR BRACKET - rear	11
---	HARNESS - erase ay	13
---	ADJUSTMENT TOOL - shutter AY	14

Numerical Index

Part No.	Description	Figure No.
---	ADJUSTMENT TOOL - clamp alignment	14
---	ADJUSTMENT TOOL - motor AY	14
---	ADJUSTMENT TOOL - sensor S3	14
---	ADJUSTMENT TOOL - bar code label	14
---	ADJUSTMENT TOOL - fold mirror	14
---	ADJUSTMENT TOOL - static brush	14
---	STAND AY - monitor	15
---	BIN - cassette storage	15
---	SHELF - monitor	15
---	LEVELING FOOT - 1.50 x 2.24 in.	15
---	BASE - floor stand	15
---	BRACKET - floor stand	15
---	MONITOR	16
---	COMPUTER BRACKET	16
---	FRONT COMPUTER BRACKET	16
---	THUMBSCREW - computer bracket	16
---	MOUNTING BRACKET - wall stand	16
---	WASHER	16
---	SCREW	16
---	WALL STAND BRACKET	16
---	BRACKET - bar code reader	16
---	KEYBOARD - standard	17
1C0554	TRANSFORMER - erase ay	8
1C7904	LATCH AY - cassette	12
1F2412	BAR CODE READER - internal	1
1F2423	LABEL - 35 x 43 cm, LLI, 136 speed	12
1F2424	LABEL - 15 x 30 cm, dental, 136 speed	12
1F2450	LABEL - 24 x 30 cm, HR-M, 68 speed	12
1F2451	LABEL - 18 x 24 cm, HR-M, 68 speed	12
1F4044	FOAM - pmt	4
1F5679	FILTER - laser, ND laser safety	2
1F5817	PLUSH - 750, left	11
1F6158	PLUSH - light seal, front	10

Numerical Index

Part No.	Description	Figure No.
1F6393	SPRING - extension, anti-reverse ay	9
1F7768	WEAR BAR - clamp, right	11
1F8430	GALVO AY	2
324259	SCREW - <i>Phillips</i> , M3 x 0.5 x 5 mm	8
365159	MOUNT - wire tie	8
3E1274	GASKET - thermal	1
3E1275	ROTOR AY	1
3E1359	RELAY	3
3H3775	LATCH AY - erase	6
477506	CARRIER - fuse	3
597341	FUSE - 8A, 250 V, slow blow	3
5E1927	COUPLING - rotor	1
5H0495	LABEL - 35 x 35 cm, GP+, 136 speed	12
5H0496	LABEL - 35 x 43 cm, GP+, 136 speed	12
5H7918	PMT AY	4
699577	SENSOR - detector	10
699579	SENSOR - emitter	10
6E5215	LIGHT SEAL - PLUSH, rear	11
6E5893	LABEL - tube side	12
6H1776	FILTER - fan	4
731068	WASHER - m3	8
7C6676	LABEL - 18 x 24 cm, GP, 136 speed	12
7C6677	LABEL - 24 x 30 cm, GP, 136 speed	12
7C6678	LABEL - 35 x 43 cm, GP, 136 speed	12
7C6679	LABEL - 18 x 24 cm, HR, 68 speed	12
7C6680	LABEL - 24 x 30 cm, HR, 68 speed	12
7E2130	LABEL - storage phosphor warning	12
7E9471	ERASE LAMP	8
851123	SCREW - 10 - 32 x 0.375 in.	10
852961	LOCK WASHER - light	10
914047	LOCK NUT - hex	8, 15
943134	SCREW - Torx m4 x 8 mm	8
943366	HOLDER - fuse	3

Numerical Index

Part No.	Description	Figure No.
955815	POWER CORD - US&C hospital grade	17
955887	SPRING - latch ay, cassette	12
9B8733	LABEL - 35 x 35 cm, GP, 136 speed	12
9F3239	LABEL - 18 x 24 cm, EHR, 80 speed	12
9F3240	LABEL - 24 x 30 cm, EHR-M, 80 speed	12
9F4123	BOARD - slow scan, coil	1
K1J7774	SENSOR - GAP, wide	11
SP182590	WIRE TIE	8
SP325608	DRIVE BELT - cassette handling ay	10
SP3B1966	POWER CORD - europe	17
SP4B0826	POWER CORD - UK	17
SP5E5632	POWER CORD - japan	17
SP5E6714	POWER CORD - china	17
SP5F7653	SUPPORT - lamp	8
SP6H3359	EXTRACTION BAR AY	1
SP6H3377	BOARD - mini mcb	3
SP6H4562	PLATE - slider, extraction bar	1
SP6H4565	PLATE - hook, extraction bar	1
SP6H4575	MOTOR AY - y-hook	1
SP6H4590	COVER - right	5
SP6H4591	COVER - left	5
SP6H4656	SPRING - compression, extraction bar	1
SP6H4957	COVER - front	5
SP6H4975	EMI PANEL AY - front	6
SP6H5190	LIGHT SHUTTER AY - right	11
SP6H5260	PLUSH AY - left	10
SP6H5471	COVER - top	5
SP6H5474	NOSE PIECE - top cover	5
SP6H5836	SHUTTER AY - front	10
SP6H6063	DRIVE AY	10
SP6H6423	SPRING - torsion, right plush	11
SP6H6527	CASSETTE GUIDE - right	5
SP6H6528	CASSETTE GUIDE - left	5

Numerical Index

Part No.	Description	Figure No.
SP6H6544	RAMP - entrance	10
SP6H6768	CAP - extraction bar	1
SP6H6814	REAR CLIP AY	11
SP6H6991	WEAR BAR - clamp, left	11
SP6H7815	SPRING - extension, extraction bar	1
SP6H7880	LCD MONITOR - 19 in., touchscreen	17
SP6H7883	LCD MONITOR - 17 in.	17
SP6H7946	ANTI-REVERSE AY - plate support	9
SP6H8005	GASKET - EMI, erase, long	6
SP6H8007	GASKET - EMI, erase, short	6
SP6H8029	CABLE - galvo motor	13
SP6H8030	CABLE - PMT	13
SP6H8080	SWITCH - main power	3
SP6H8426	BAFFLE - afterglow	4
SP6H8543	CABLE - laser control	13
SP6H8593	CABLE - ethernet	13
SP6H8594	CABLE - ieb board to mini-mcb board	13
SP6H9388	STATIC BRUSH - extraction bar	1
SP6H9581	COMPUTER - model 8810	17
SP6H9808	EXTENSION SPRING - clutch	9
SP6H9897	FAN AY - collector	4
SP6H9899	SENSOR AY - plate present	1
SP748360	SCREW - <i>Torx</i> , m4 x 12 mm	1. 8
SP748361	SCREW - <i>Sems</i> , hex flange head	2
SP7E7744	SCREW - <i>Sems</i> , hex, m4 x 10 mm	4
SP7H2391	DRIVE MOTOR AY - cassette handling ay	10
SP7H2639	POWER SUPPLY - 19 in. monitor	17
SP7H2640	CABLE - VGA	13
SP7H2641	CABLE - 19 in. monitor, USB	17
SP7H2642	CABLE - 19 in. monitor, audio	17
SP7H2685	ERASE AY - complete	8
SP7H2777	HARNESS - power supply	3
SP7H2860	SUPPORT - transformer	3

Numerical Index

Part No.	Description	Figure No.
SP7H2956	BOOT - bar code reader, wall stand	16
SP7H2968	BAFFLE AY - erase ay	8
SP7H3043	TRANSFORMER - isolation	3
SP7H3310	POWER SUPPLY	3
SP7H3363	CASSETTE HANDLING AY - complete	10
SP7H3840	SPRING - light seal ay, left	10
SP7H3855	BOARD - PMT, 1.0 mz	4
SP7H3856	CONTROL DISPLAY AY	5
SP7H3964	STAND - 19 in. monitor	17
SP7H3991	PLUSH - 750, right	11
SP7H4100	PLUSH - clamp bar	11
SP7H4386	SCREW - captive	5
SP7H4840	HARNESS - mini-mcb to cassette handling ay	10. 13
SP7H4879	TRAY - keyboard, floor stand	15
SP7H5083	TRAY - keyboard, wall stand	16
SP7H5087	KEYBOARD - mini	17
SP7H5089	ARM - bar code	16
SP7H5101	LINER AY - bin, floor stand	15
SP7H5107	MONITOR ARM - wall stand	16
SP7H5399	BAR CODE READER - handheld	16. 17
SP7H5400	USB CABLE - bar code reader	16. 17
SP7H5401	STAND - handsfree bar code reader	17
SP7H5403	CUP - bar code reader, desktop	17
SP7H5420	FOLD MIRROR AY	4
SP7H5533	BOARD - IEB	3
SP7H5992	LIGHT SEAL AY - right	5
SP7H5993	LIGHT SEAL AY - left	5
SP7H6000	CIRCUIT BREAKER - power supply	3
SP7H6068	SLIDE AY - slow scan, complete	1
SP7H6142	COLLECTOR AY - complete	4
SP7H6238	MOTOR AY - extraction bar	1
SP7H6652	SUPPORT AY - plate	7
SP7H6765	SPRING - extension, clamp	10

03DEC07
IP5258-1
Page
44 of 53

Numerical Index

Part No.	Description	Figure No.
SP7H7434	MINI FAST SCAN	2
SP7H7435	MINI FAST SCAN - complete	2
SP943600	SCREW - <i>Torx</i> , m4 x 18 mm	1. 4
SP9B9445	SCREW - <i>Torx</i> m2.5 x 10 mm	8
SP9F0341	MOTOR - cassette handling ay, DC	10
SPSP2B4081	SCREW - <i>Torx</i> , m4 x 8 mm	5
TL3E1279	SWITCH - optical	10
TL7H4608	ADJUSTMENT TOOL KIT	14

Section 3: Alphabetical Index

Alphabetical Index

Part No.	Description	Figure No.
---	ADJUSTMENT TOOL - bar code label	14
---	ADJUSTMENT TOOL - clamp alignment	14
---	ADJUSTMENT TOOL - fold mirror	14
---	ADJUSTMENT TOOL - motor AY	14
---	ADJUSTMENT TOOL - sensor S3	14
---	ADJUSTMENT TOOL - shutter AY	14
---	ADJUSTMENT TOOL - static brush	14
TL7H4608	ADJUSTMENT TOOL KIT	14
SP6H7946	ANTI-REVERSE AY - plate support	9
SP7H5089	ARM - bar code	16
---	ARM - plate support	7
---	ARM BOOSTER - plate support	9
SP6H8426	BAFFLE - afterglow	4
SP7H2968	BAFFLE AY - erase ay	8
---	BALL BEARING - flanged, 0.25 x 0.37 in.	7
SP7H5399	BAR CODE READER - handheld	16. 17
1F2412	BAR CODE READER - internal	1
---	BASE - floor stand	15
---	BEARING - needle	1
---	BIN - cassette storage	15
SP7H5533	BOARD - IEB	3
SP6H3377	BOARD - mini mcb	3
SP7H3855	BOARD - PMT, 1.0 mz	4
9F4123	BOARD - slow scan, coil	1
---	BOLT - shoulder	1
SP7H2956	BOOT - bar code reader, wall stand	16
---	BRACKET - bar code reader	16
---	BRACKET - disc, friction	9
---	BRACKET - floor stand	15
---	BRACKET - plush, rear	11
---	BRACKET - power supply	3

Alphabetical Index

Part No.	Description	Figure No.
SP7H2642	CABLE - 19 in. monitor, audio	17
SP7H2641	CABLE - 19 in. monitor, USB	17
SP6H8593	CABLE - ethernet	13
SP6H8029	CABLE - galvo motor	13
SP6H8594	CABLE - ieb board to mini-mcb board	13
SP6H8543	CABLE - laser control	13
SP6H8030	CABLE - PMT	13
SP7H2640	CABLE - VGA	13
SP6H6768	CAP - extraction bar	1
477506	CARRIER - fuse	3
SP6H6528	CASSETTE GUIDE - left	5
SP6H6527	CASSETTE GUIDE - right	5
SP7H3363	CASSETTE HANDLING AY - complete	10
SP7H6000	CIRCUIT BREAKER - power supply	3
---	CLUTCH - 22 x 8 x 9 mm	9
---	COLLECTOR AY - 2 PMT	4
SP7H6142	COLLECTOR AY - complete	4
---	COMPRESSION DISC - outer	9
---	COMPRESSION DISK - inner	9
SP6H9581	COMPUTER - model 8810	17
---	COMPUTER BRACKET	16
---	CONNECTOR	3
SP7H3856	CONTROL DISPLAY AY	5
5E1927	COUPLING - rotor	1
SP6H4957	COVER - front	5
SP6H4591	COVER - left	5
---	COVER - optics ay	2
---	COVER - pmt	4
---	COVER - relay	3
SP6H4590	COVER - right	5
SP6H5471	COVER - top	5
SP7H5403	CUP - bar code reader, desktop	17
---	DISC AY - plate support	9

Alphabetical Index

Part No.	Description	Figure No.
SP6H6063	DRIVE AY	10
SP325608	DRIVE BELT - cassette handling ay	10
SP7H2391	DRIVE MOTOR AY - cassette handling ay	10
---	EMI PANEL - left	6
---	EMI PANEL - rear	6
---	EMI PANEL - right	6
SP6H4975	EMI PANEL AY - front	6
SP7H2685	ERASE AY - complete	8
7E9471	ERASE LAMP	8
---	EXTENSION SPRING - arm booster	9
SP6H9808	EXTENSION SPRING - clutch	9
SP6H3359	EXTRACTION BAR AY	1
SP6H9897	FAN AY - collector	4
6H1776	FILTER - fan	4
1F5679	FILTER - laser, ND laser safety	2
---	FILTER BRACKET - laser	2
---	FLAT WASHER - 0.375 x 0.170 x 0.024 in.	1
1F4044	FOAM - pmt	4
---	FOLD MIRROR	4
SP7H5420	FOLD MIRROR AY	4
---	FRAME - erase ay	8
---	FRONT COMPUTER BRACKET	16
597341	FUSE - 8A, 250 V, slow blow	3
1F8430	GALVO AY	2
SP6H8005	GASKET - EMI, erase, long	6
SP6H8007	GASKET - EMI, erase, short	6
3E1274	GASKET - thermal	1
---	GROUND STRAP	1
---	HARNESS - ballasts and sockets	8
---	HARNESS - erase ay	13
SP7H4840	HARNESS - mini-mcb to cassette handling ay	10, 13
---	HARNESS - plug and sensor	8
SP7H2777	HARNESS - power supply	3

Alphabetical Index

Part No.	Description	Figure No.
943366	HOLDER - fuse	3
---	HOUSING - PMT	4
SP7H5087	KEYBOARD - mini	17
---	KEYBOARD - standard	17
1F2424	LABEL - 15 x 30 cm, dental, 136 speed	12
9F3239	LABEL - 18 x 24 cm, EHR, 80 speed	12
7C6676	LABEL - 18 x 24 cm, GP, 136 speed	12
7C6679	LABEL - 18 x 24 cm, HR, 68 speed	12
1F2451	LABEL - 18 x 24 cm, HR-M, 68 speed	12
9F3240	LABEL - 24 x 30 cm, EHR-M, 80 speed	12
7C6677	LABEL - 24 x 30 cm, GP, 136 speed	12
7C6680	LABEL - 24 x 30 cm, HR, 68 speed	12
1F2450	LABEL - 24 x 30 cm, HR-M, 68 speed	12
5H0495	LABEL - 35 x 35 cm, GP+, 136 speed	12
9B8733	LABEL - 35 x 35 cm, GP, 136 speed	12
5H0496	LABEL - 35 x 43 cm, GP+, 136 speed	12
7C6678	LABEL - 35 x 43 cm, GP, 136 speed	12
1F2423	LABEL - 35 x 43 cm, LLI, 136 speed	12
---	LABEL - erase ay	8
7E2130	LABEL - storage phosphor warning	12
6E5893	LABEL - tube side	12
1C7904	LATCH AY - cassette	12
3H3775	LATCH AY - erase	6
SP6H7883	LCD MONITOR - 17 in.	17
SP6H7880	LCD MONITOR - 19 in., touchscreen	17
---	LEAF SPRING - anti-reverse	9
---	LEVELING FOOT - 1.50 x 2.24 in.	15
6E5215	LIGHT SEAL - PLUSH, rear	11
SP7H5993	LIGHT SEAL AY - left	5
SP7H5992	LIGHT SEAL AY - right	5
SP6H5190	LIGHT SHUTTER AY - right	11
SP7H5101	LINER AY - bin, floor stand	15
914047	LOCK NUT - hex	8. 15

Alphabetical Index

Part No.	Description	Figure No.
852961	LOCK WASHER - light	10
SP7H7434	MINI FAST SCAN	2
SP7H7435	MINI FAST SCAN - complete	2
---	MONITOR	16
SP7H5107	MONITOR ARM - wall stand	16
SP9F0341	MOTOR - cassette handling ay, DC	10
SP7H6238	MOTOR AY - extraction bar	1
SP6H4575	MOTOR AY - y-hook	1
---	MOTOR GUARD - slow scan	1
---	MOUNT - collector, front	4
---	MOUNT - collector, rear	4
---	MOUNT - flat cable	4
---	MOUNT - fold mirror	4
365159	MOUNT - wire tie	8
---	MOUNTING BRACKET - wall stand	16
SP6H5474	NOSE PIECE - top cover	5
---	PIVOT BRACKET - plate support	7
SP6H4565	PLATE - hook, extraction bar	1
SP6H4562	PLATE - slider, extraction bar	1
---	PLUG	3
1F5817	PLUSH - 750, left	11
SP7H3991	PLUSH - 750, right	11
SP7H4100	PLUSH - clamp bar	11
1F6158	PLUSH - light seal, front	10
SP6H5260	PLUSH AY - left	10
5H7918	PMT AY	4
---	POST PIVOT - eccentric	9
SP5E6714	POWER CORD - china	17
SP3B1966	POWER CORD - europe	17
SP5E5632	POWER CORD - japan	17
SP4B0826	POWER CORD - UK	17
955815	POWER CORD - US&C hospital grade	17
SP7H3310	POWER SUPPLY	3

Alphabetical Index

Part No.	Description	Figure No.
SP7H2639	POWER SUPPLY - 19 in. monitor	17
---	PULLEY - drive motor, 20T	10
---	PUSHER ARM - plate support	9
---	RACK AY - power supply	3
SP6H6544	RAMP - entrance	10
SP6H6814	REAR CLIP AY	11
3E1359	RELAY	3
---	RIVET - 0.123 x 7/32 in.	7
3E1275	ROTOR AY	1
---	SCREW	16
851123	SCREW - 10 - 32 x 0.375 in.	10
SP7H4386	SCREW - captive	5
324259	SCREW - <i>Phillips</i> , M3 x 0.5 x 5 mm	8
SP748361	SCREW - <i>Sems</i> , hex flange head	2
SP7E7744	SCREW - <i>Sems</i> , hex, m4 x 10 mm	4
SP9B9445	SCREW - <i>Torx</i> m2.5 x 10 mm	8
943134	SCREW - <i>Torx</i> m4 x 8 mm	8
SP748360	SCREW - <i>Torx</i> , m4 x 12 mm	1. 8
---	SCREW - <i>Torx</i> , m4 x 18 mm	6
SP943600	SCREW - <i>Torx</i> , m4 x 18 mm	1. 4
---	SCREW - <i>Torx</i> , m4 x 6 mm	9
SPSP2B4081	SCREW - <i>Torx</i> , m4 x 8 mm	5
---	SEAL - light-power switch	6
699577	SENSOR - detector	10
699579	SENSOR - emitter	10
K1J7774	SENSOR - GAP, wide	11
SP6H9899	SENSOR AY - plate present	1
---	SENSOR BRACKET - rear	11
---	SHELF - monitor	15
---	SHOULDER SCREW - 10-32 x 0.25 x 1.00 in.	7
---	SHOULDER SCREW - M4	9
---	SHOULDER SCREW - M4 x 5.0 x 25 mm	9
SP6H5836	SHUTTER AY - front	10

Alphabetical Index

Part No.	Description	Figure No.
---	SKID - arm, plate support	7
---	SLIDE - encoder AY	1
---	SLIDE AY - slow scan	1
SP7H6068	SLIDE AY - slow scan, complete	1
SP6H4656	SPRING - compression, extraction bar	1
1F6393	SPRING - extension, anti-reverse ay	9
SP7H6765	SPRING - extension, clamp	10
SP6H7815	SPRING - extension, extraction bar	1
955887	SPRING - latch ay, cassette	12
SP7H3840	SPRING - light seal ay, left	10
SP6H6423	SPRING - torsion, right plush	11
---	SPRING PIN - 0.094 x 0.812 in., steel	10
SP7H3964	STAND - 19 in. monitor	17
SP7H5401	STAND - handsfree bar code reader	17
---	STAND AY - monitor	15
---	STANDOFF	3
---	STANDOFF - clinch, 10-32	7
SP6H9388	STATIC BRUSH - extraction bar	1
SP5F7653	SUPPORT - lamp	8
SP7H2860	SUPPORT - transformer	3
SP7H6652	SUPPORT AY - plate	7
SP6H8080	SWITCH - main power	3
TL3E1279	SWITCH - optical	10
---	THUMBSCREW - computer bracket	16
---	THUMBSCREW - m3	2
1C0554	TRANSFORMER - erase ay	8
SP7H3043	TRANSFORMER - isolation	3
SP7H4879	TRAY - keyboard, floor stand	15
SP7H5083	TRAY - keyboard, wall stand	16
SP7H5400	USB CABLE - bar code reader	16, 17
---	WALL STAND BRACKET	16
---	WASHER	16
---	WASHER - 0.62 x 0.255 x 0.062 in.	7, 9

03DEC07
IP5258-1
Page
52 of 53

Alphabetical Index

Part No.	Description	Figure No.
---	WASHER - curved, 0.49 x 0.25 x 0.011 in.	7
---	WASHER - lock	8
---	WASHER - lock, m3, internal tooth	8
731068	WASHER - m3	8
SP6H6991	WEAR BAR - clamp, left	11
1F7768	WEAR BAR - clamp, right	11
SP182590	WIRE TIE	8

Section 4: Publication History

Publication Date	Publication No.	ECO No.	Changed Pages	File Name	Description
26NOV07	IP5258-1	---	---	ip5258_1.fm	New Publication
03DEC07	IP5258-1	---	6	ip5258_1.fm	Revision

DirectView is a trademark of Carestream Health.
Kodak is a trademark of Kodak used under license.